

SemperSure[®]

NITRILE EXAMINATION GLOVES

SuperCare[™]
for your Hands



Skin-Friendly Examination Nitrile Gloves

- **Thickness:** Competes with a 3-4 mil glove
- Manufactured without accelerators
- Non-chlorinated
- Low dermatitis potential
- Reduced risk of skin irritations
- Excellent barrier and chemical resistance
- Tested for use against Fentanyl
- Meets USP 800 requirements for PPE compliance for chemotherapy
- Tested for protection against gastric acid
- Meets FDA 21 CFR Part 177 requirements for food handling

RECOMMENDED USES



PRODUCT SPECS



200 / 10 /

(180 XL)

REORDER INFO		Glove Color: Grape	
SUNF201	X-Small	200 gloves/box	10 boxes/case
SUNF202	Small	200 gloves/box	10 boxes/case
SUNF203	Medium	200 gloves/box	10 boxes/case
SUNF204	Large	200 gloves/box	10 boxes/case
SUNF205	X-Large	180 gloves/box	10 boxes/case

SAMPLES: Available by request - ask your sales representative

For maximum protection, we recommend changing nitrile gloves every 20-30 minutes

For more information about protection against fentanyl, please contact your sales rep.

Rev01202023

